
D-Heart Portable ECG Device

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CONTENT

English	4
Italiano	30
Français	58
Deutsch	84
Español	112
Português	140

EN

D-HEART PORTABLE ECG DEVICE USER MANUAL

1.	CONTENT OF THE PACKAGE AND SERVICE LIFE	6	7.	HARDWARE PLATFORM AND OPERATING SYSTEM REQUIREMENTS	15
2.	IMPORTANT WARNINGS	6	8.	ADDITIONAL ACCESSORY CLOUD	16
3.	INDICATIONS FOR USE	7	9.	INTERPRETATION OF RESULTS	17
4.	DEVICE OPERATION	8	10.	CLINICAL PERFORMANCE TESTING	17
	a. POWER ON/OFF	8	11.	ENVIRONMENTAL CONDITIONS	17
	b. DEVICE STATUS	8	12.	DECLARATION OF CONFORMITY	18
5.	HOW TO USE THE D-HEART DEVICE AND APP	9	13.	INTENDED USE, INTENDED USER, CONTRA-INDICATIONS AND PRECAUTIONS, RESIDUAL RISKS	18
	5.1 APP DOWNLOAD	9	14.	MAINTENANCE, CLEANING AND TROUBLESHOOTING	21
	5.2 REGISTRATION	9		a. CLEANING AND STERILIZING	21
	5.3 AUTHENTICATION	9		b. MAINTENANCE	21
	5.4 SETTINGS	9		c. TROUBLESHOOTING	21
	5.5 BLUETOOTH PAIRING	10	15.	LEGENDA OF SYMBOLS	23
	5.6 HOW TO REGISTER AN ECG	10	16.	SUPPORT IN CASE OF NEED AND ADVERSE EVENT REPORTING	25
	5.7 ELECTRODE POSITIONING HELP	10	17.	ELECTROMAGNETIC COMPATIBILITY	25
	5.8 RECORDING THE ECG	12	18.	DISPOSAL OF THE DEVICE	28
	5.9 SAVING THE ECG	12	19.	CONTACT US	29
	5.10 RESULTS	12			
	5.11 EXAMS	14			
	5.12 PATIENTS	14			
6.	OUTPUTS	14			

ATTENTION: Read this manual carefully before using the instrument.
Please keep it in a safe place for future reference.

1. CONTENT OF THE PACKAGE AND SERVICE LIFE

- **D-Heart Portable ECG Device** with lanyard¹
- Charging cable
- Soft case with embedded Wireless charger
- Quick Guide and instructions for use
- Additional label warning users not to lose the Operator's Manual and to read Manual before use

The D-Heart device service life is estimated to be 4 years.

The service life of the charging cable and housing is estimated at 5 years.

Electrodes, not included in the package, are disposable (and have an expiry date, indicated on their packaging).

2. IMPORTANT WARNINGS

- For a correct use of the D-Heart, please read all the instructions included in this manual before using the device
- MR unsafe! Do not expose the device to a magnetic resonance (MR) environment
- The D-Heart is not defibrillator-proof. The user must remove the D-Heart device from the body before using any high voltage device such as an AED (Automated External Defibrillator)
- The device is not to be used on people with pacemakers
- Do not use the device while charging it – The device will lose Bluetooth connection with the smartphone while charging
- The device is intended for adults. Keep the device out of reach of children, up high, out of sight, and out of reach to prevent potential ingestion hazards of small parts (e.g. disposable electrodes for use with the device)
- The D-Heart App intended use is not for direct diagnosis or for monitoring vital physiological processes where the nature of variations is such that it could result in immediate danger to the patient
- User authentication is mandatory to access the App and consequently use the device. Labeling advises users on the correct password creation and on the correct password update and management
- For Health Professionals: there may be patients of the same name, be careful
- Labeling advises users on the correct management of their own

¹ In case of need, the holder of the lanyard can be detached from the D-Heart Portable ECG Device.

- smartphones and device to protect them from theft
- Protect the web application from hacker attacks with firewalls and antivirus and any other state-of-the-art steps
- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the D-Heart, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could occur
- The device shall not be operated in combination with high frequency surgical equipment
- The device shall not be operated in combination with flammables
- The device shall not be operated in combination with flammable anesthetics
- The device shall not be operated in environments saturated by oxygen
- The storage range is narrower than the range indicated in IEC 60601-1-11+Ed.2.0(2015)+Amd1.(2020). Make sure to storage between -20°C and +30°C
- Trace acquisition via the D-Heart device can be performed by any user, autonomously in both home and hospital environments. This is possible thanks to the electrode positioning guide on the D-Heart App and the design of the device that allows the cables to be correctly unwound without the risk of reversing them.

3. INDICATIONS FOR USE

D-Heart is a portable ECG. The device is designed to record, store and transfer the rhythms of an eight- and twelve-channel electrocardiogram (ECG). In the eight-channel mode, D-Heart can record DI, DII, V2 and V5 unipolar simultaneously, and derive aVR, aVF and aVL. In the twelve-channel mode, D-Heart can simultaneously record DI, DII, V2, and V5 unipolar simultaneously, derive aVR, aVF, and aVL, and non-simultaneously record V1, V3, V4, and V6 unipolar. The device is intended to be operated by healthcare professionals and patients with known or suspected heart conditions. The evaluation and reporting of electrocardiographic tracing, on the other hand, should only and exclusively be carried out by qualified health professionals. The device can be

operated in hospitals, general practitioner offices, pharmacies, and in home locations. The device is intended for adults. The D-Heart device cannot be used in ambulatory conditions; the device is not an ambulatory monitor as the user/patient needs to stay still when registering the exam. The D-Heart is a portable ECG device, meaning that it can be carried and used everywhere (e.g., carried inside a luggage when travelling, inside a bag or backpack when walking to be used in different home environments, office environments, outside environments, etc.), however, the user/patient cannot move when registering the ECG.

3.1 INTENDED PURPOSE

The D-Heart device is a portable ECG device consisting of a Bluetooth front-end and an APP running on Android and iOS, designed to help patients to detect pathological states of the heart through the presentation of the heart's electrical signal and some heart parameters (BPM, PR, QRS intervals).

The D-Heart App allows to:

- Display 8 or 12 leads ECGs in real time and some heart parameters (BPM, PR, QRS intervals) through the interpretation of data acquired by the Bluetooth front end.
- Send the file with the ECG display via https for use in other ECG analysis software.

The D-Heart device is intended for use on people over 18 years of age.

The D-Heart device is intended to be used by health professionals and patients in both domestic and hospital environments.

4. DEVICE OPERATION

a. POWER ON/OFF

- When the device is OFF (Led off), press the central button to switch the device ON.
- When the device is ON (Led blinking), press the central button to switch the device OFF.

b. DEVICE STATUS

- OFF (Led off)
- ON/NOT PAIRED
Led single blinking; timing: 1s
- ON/CONNECTED
Led ON (fixed)
- ON/MEASURING
Led double blinking; timing: 1s
- ON/BATTERY DISCHARGED
Led blinking 3 times; then OFF
- ON/FIRMWARE UPDATE
Led fast blinking; timing: 0,3s
- OFF/CHARGING
Led single blinking; timing: 5s
- OFF/CHARGED
Led single blinking; timing: 1,5s

5. HOW TO USE THE D-HEART DEVICE AND APP

5.1 APP DOWNLOAD

The device can only be used together with the D-Heart App. The D-Heart App can be downloaded free of charge from Google Play Store when using Android devices and the App Store when using iOS-Apple devices.

5.2 REGISTRATION

Once the App is downloaded the user will need to register and create a new profile.

The password must be at least 8 characters long, not more than 20, with an uppercase character, a lower case character, a number and a special character among the following (@!\$%*?&).

Passwords must be kept safe and must not be shared with anyone or any other third party. Consider changing your password cyclically (e.g. every 6 months).

5.3 AUTHENTICATION

Only authorized users can access the App. The user can authenticate, enable Touch ID login or create a new account (after acknowledging Terms and Conditions). After the authentication, if the account is not yet completed, the application will prompt the user to complete their profile.

5.4 SETTINGS

The settings page (top right corner of the homepage) contains specific application settings, the opportunity to choose the App language, ECG settings, Device and Bluetooth settings, the opportunity to export the ECGs already registered, FAQs, and last accepted Terms and Conditions.

The Terms and Conditions section also contains the App Label.

The Device and Bluetooth Settings section allows to see the device firmware and hardware versions and the Serial Number of the device. In order to obtain this information, it is necessary to connect the device to the App (please see 5.5 Bluetooth pairing on how to connect the device to the App). The FAQs section contains all the instructions to correctly record an ECG.

5.5 BLUETOOTH PAIRING

This function guides the user to the connection of the Smartphone to the D-Heart ECG device.

In order to connect the device to the smartphone, click on the device icon on the top left side of the App homepage and then click 'Begin'. Click on DHRT and the device will be connected to the smartphone. Alternatively, click on the settings page, then click on Device and Bluetooth settings and then click on Start Searching. Click on DHRT and the device will connect to the smartphone. After the first connection, the saved device will be automatically connected every time the application starts. If the device search does not find any D-Heart available, the App shows a button to contact the D-Heart support team.

5.6 HOW TO REGISTER AN ECG

The D-Heart ECG is NOT a cardiac holter. In order to register the ECG, the user needs to be sitting down and to stay still. The App will communicate to the user if they are moving too much. Before registering the ECG, make sure not to be close to other electronic equipment such as TV screens, PCs, radios, electric beds, etc. Which may cause interference in the Bluetooth signal. In order to register an ECG, click on the orange heart button in the center bottom of the screen. The App will ask how many leads the user wants to record. Select either 8 or 12 leads. The app will automatically show a tutorial video (one for the 8 lead and a different one for the 12 lead) to first time users. The videos explain and show the procedure to correctly register an ECG with the device. The user can choose from the App whether to acquire an 8- or 12-lead trace. While the 8-lead trace is easy to acquire by lay users independently under the guidance of the D-Heart App, it is recommended that the 12-lead trace should be acquired by experienced users only or by patients supported by a caregiver or healthcare professional.

5.7 ELECTRODE POSITIONING HELP

Follow the instructions below carefully: Incorrect positioning of the electrodes may result in the ECG waveform not being recorded or in unreadable or incorrect signals being recorded. In order to guide the user to the correct positioning of the electrodes on the body to record an 8-lead ECG by non-professional users, the camera of the smartphone uses a powerful electrode positioning feature. The user is invited to take a picture of their chest and the augmented reality feature automatically produces a picture showing the correct position of the electrodes (see Fig. 9). In order to take the picture, the user needs to be sitting down and to fit their bare torso inside the grey shape showed by the App. There is no fixed distance to sit away from the camera

to take the placement photo, but the user needs to see themselves inside the grey shape. After connecting the electrodes, following the guidelines shown in the resulting photo, you can start with the recording.


Fig. 9: the electrode placement assistance flow

If a 12-lead ECG is selected by the user, the same tutorial ECG will be opened and standard D-Heart ECG device positioning will be suggested. To record the other precordial leads, 4 sequential animated GIFs will guide the user on the correct precordial electrode placement as described in Fig. 10.


Fig. 10: electrode positioning assistance flow for 12-lead ECGs

5.8 RECORDING THE ECG

During the registration, the user can see the 8-lead ECG on the screen (Lead DI, DII, DIII, aVR, aVL, aVF, V2 and V5, Fig. 2). There will be a countdown that shows the remaining registration time and the real time Beat per Minute (BPM). If there are errors in the communication (e.g. Packet loss with the Bluetooth device), the registration will restart automatically. If the D-Heart ECG disconnects during the exam, the exam restarts from the electrode positioning help.

To record the 12-lead ECG, the user will be asked to record a standard 8-lead ECG as described above. Subsequently, following the first animated GIF, the App will ask to move the black precordial electrode to position V6. The ECG will restart displaying the 8-lead ECG (same Lead DI, DII, DIII, aVR, aVL, aVF, V2 and now V6) and record only the V6 lead tracing. The process will then be repeated for all the remaining precordial Leads (V4, V3, V1), as shown by Fig. 11. A final PDF showing the 12-lead ECG will be displayed on the screen.

Device operation is also described in the *Quick Guide for D-Heart Portable ECG Device* included in the package and in the Video Tutorials available at www.d-heartcare.com.

5.9 SAVING THE ECG

Once the ECG has been registered, the user will need to associate it to a patient by clicking on 'No Patient Selected' and select a pre-existing patient or create a new one. The user will then need to select the symptom by clicking on 'Symptom'. The user will also have the opportunity to insert an additional note by clicking on 'Note'. The user will then be able to save the ECG by clicking on 'Save'.

5.10 RESULTS

Once the ECG is saved the user should scroll down the page of the exam page and will be able to download or directly share or print the registered ECG by clicking on 'PDF'. Laymen will not be able to read and analyze the ECG, but they can choose to share it with a clinician, their General Practitioner or their cardiologist to receive an analysis of the ECG, thus a result. The acquisition of the trace takes place in total autonomy without the intervention of a professional. The latter may have to be involved at a later stage for reading and interpreting the ECG. In order to receive an analysis and interpretation of the ECG from D-Heart, the user can also choose to use the additional tele-cardiology


Fig. 11: 12 Lead ECG Recording Flowchart

service by clicking on 'Send to Telecardiology'.

The user will be shown the number of credits remaining to send the ECG to Telecardiology under the 'Send to Telecardiology' button. Additional credits can be bought on the D-Heart e-commerce store <https://www.d-heartcare.com/single-telecardiology-report-credit.html>

If the user chooses to send the ECG to the telecardiology, they will receive a report signed by a cardiologist back within 15 minutes.

The telecardiology service is an additional service and is not mandatory. People who do not have the necessary medical knowledge and who

choose not to share their ECG with their cardiologist, nor to use the additional telecardiology service, will not be able to receive an interpretation of their exam.

5.11 EXAMS

This function shows the ECG archive or the list of the last exams registered (Exams icon on the right bottom side of the homepage). Every exam will be shown with the date, the exam results and the patient's name. Filters are available to select a specific exam from the list.

5.12 PATIENTS

The patients' section (Patients icon on the left bottom side of the homepage) allows users to see all patients' profiles saved in the App and allows to add new patients' profiles. Patients will be displayed in alphabetical order.

If there are more than one patient of the same name, make sure to open the patient information sheet to check that the correct patient is selected. Pay special attention to the patient's date of birth as a way to distinguish homonymous patients.

6. OUTPUTS

Current device outputs are considered non-obvious measurements. The subject device provides 8-Lead ECGs of 12, 24, 48 or 60 seconds and 12-lead ECGs of 6 seconds each. It is possible to change the track acquisition time, directly from the application.

The subject device also automatically measures the intervals PR, QRS, and BPM:

- The interval PR is defined as the time calculated in milliseconds that the cardiac impulse takes to travel from its beginning (in the sino atrial node) up to the atrioventricular node.
- The QRS interval measures the time calculated in milliseconds that the cardiac impulse takes to depolarize the left ventricle.
- The BPM refers to cardiac frequency during 60 seconds.

The ECG output should be interpreted by a health professional specifically trained in ECG interpretation.

The interpretation can be performed directly on the device, remotely

by sharing the ECG via a third party telecardiology platform or it can be printed if the smartphone can be connected to a printer via Bluetooth.

7. HARDWARE PLATFORM AND OPERATING SYSTEM REQUIREMENTS

The D-Heart App requires at least 2 GB of memory, a BLE version 4.0, a display of minimum size 5" and resolution 720 on short screen side. The D-Heart App is tested and guaranteed to work on the following smartphones:

1. iPhone 7 and newer, with iOS version 12 and later

- 1.3 GHz Dual Core processor;
- RAM 1 Gb;
- 8 Mega Pixel rear camera;
- 3264 x 2448 ppi screen resolution;
- Front camera 1.2 mega pixels;
- Video 30 fps.

2. Android 8 and later, with a CPU Exynos 9 Octa or equivalent. Less powerful CPU are not guaranteed to work but are guaranteed to not use incomplete or malformed data.

- 1.2 GHz 4 Core processor;
- 2 GB RAM;
- 13 Mega Pixel rear camera;
- 4160 x 3120 ppi screen resolution;
- 5 Mega Pixel front camera;
- Video 30 fps.

D-Heart is compatible with various Android and iOS smartphones and tablets. For an updated list of compatible devices, please visit the D-Heart website at www.d-heartcare.com/pages/compatibilit

D-Heart is not compatible with iOS version 13.4 and all Huawei Lite smartphones. Other smartphones and tablets that are not included in the above list may work with the device, but D-Heart cannot guarantee their perfect pairing and functioning. At least on a monthly basis, check the availability of App and Firmware updates and make the required upgrades accordingly.

- Make sure to set a PIN number to unlock your phone
- Consider enabling mobile tracking services

- Consider installing an anti-theft software onto your device
- Consider using an anti-theft phone case
- Take note of your phone's unique ID number
- Overall, keep your smartphone out of sight and do not leave it unattended in public / crowded places.

8. ADDITIONAL ACCESSORY CLOUD (Not covered by certification)

The device may communicate with Clouds that are not proprietary to Strena Medical S.r.l and not covered by the product's MDR certification. If the Cloud is used, at the user's own risk, it is important that it has the following characteristics:

- Cloud must have a protocol to manage user registration and login. Password management (reset, password forgotten) must be provided.
- Cloud must have a protocol to store user data allowing login from different smartphones.
- Cloud must have a protocol to store photos acquired by mobile applications during self-positioning of electrodes.
- Cloud must have a protocol to locate and delete Users specific photos.
- Cloud must have a protocol to collect data from users and store it in a database in an anonymous form.
- Cloud must have a protocol to provide a way to connect it to a third-party cloud software for remote medical report of smartphone acquired ECG data.
- Cloud must have a protocol to manage a "credit system" for each user. "Credit" must be used to consume remote medical reports.
- Cloud must have a protocol to follow a serverless approach to maintain running costs as low as possible.
- Cloud shall be developed using AWS technologies.
- Cloud must not manage medical or sensible data.
- Cloud must have a protocol to be integrated with an external e-Commerce site, and shall allow this site to increment "Credits" for a specified user.
- Cloud must have a protocol to manage at least 1000 simultaneous connection from apps
- Cloud must have a protocol to store at least 1000 exams per user.

9. INTERPRETATION OF RESULTS

The present ECG device allows the acquisition of an 8 or 12-lead ECG.

Once registered independently, the ECG, should be interpreted by trained health professionals, persons legally authorized by the national laws and regulations, with appropriate skill on muscular anatomy and functionality. All end users who are not trained to interpret ECGs may choose to use the tele-cardiology platform to receive a medical report and interpretation or should send their ECG in PDF format to their cardiologist/health professional for a correct interpretation. Information such as ECG interval length will be automatically provided by the device.

10. CLINICAL PERFORMANCE TESTING

D-Heart has been clinically tested to prove its efficiency and quality recording for both 8 and 12 lead ECGs, even in the case of autonomous trace acquisition by lay users.

Consecutive patients with a diagnosis of Hypertrophic Cardiomyopathy referred for outpatient control at a referral national institution for cardiomyopathies were enrolled for clinical testing from April to May 2017. All patients were older than 18 years old.

D-Heart® and 12-lead ECG recordings were obtained (within 2-5 min) in each subject. ESAOTE 4230 (FDA cleared 2/2/1998) model was used as a standard ambulatory ECG device, intended as the 'gold standard'.

D-Heart® tracings proved effective and accurate in the stratification of ECG abnormalities as compared to the 12-Lead ECG. D-Heart® tracings were reliable for the identification of LVH as well as for arrhythmias such as AF and Ventricular Ectopic Beats. Accuracy of ECG intervals was comparable between D-Heart and ESAOTE 4230.

11. ENVIRONMENTAL CONDITIONS

For an optimal functioning of the device, the following environmental conditions must be followed:

- Operating Temp (battery charging and discharging): +10°C to +35°C
- Storage Temp: -20°C to +30°C
- Storage Humidity: 45% - 75% (relative)
- Atmospheric pressure range: 700 hPa to 1060 hPa

If the device is stored at a temperature between -20°C and +10°C, the user will need to wait for at least 60 minutes for the device to reach its operational temperature in the operative environment.

12. DECLARATION OF CONFORMITY

Hereby, Strena Medical S.r.l Società Benefit declares that this device, D-Heart (REF ECG000-DHR20) is in compliance with all essential requirements and other provisions laid down in EU directives: 2014/53/EU (RED), 2011/65/EU (RoHS) and Regulation (EU) 745/2017 (MDR). The full text of the EU Declaration of Conformity is available at: <https://www.d-heartcare.com/en/pages/qualita-certificazioni> In accordance with European Commission Decision N°2000/299/EC of 06/04/2000 the frequency band used by this product is harmonised in all EU countries, therefore this class I product can be freely used in all EU countries.

D-Heart Portable ECG Device is compliant with the following standards and directives/regulations:

- ETSI EN 301 489-1 V2.2.3 (2019-11)
- ETSI EN 301 489-17 V3.2.4 (2020-09)
- ETSI EN 300 328: v2.2.2 (2019-07)
- IEC 60601-1-2:2014/AMD1:2020
- IEC 60601-2-25:2011
- IEC 60601-1-11:2015+AMD1:2020
- IEC 60601-1:2005+AMD2:2020
- 2014/53/UE (RED)
- 2011/65/UE (RoHS)
- Regulation (EU) 745/2017 (MDR)

13. INTENDED USE, INTENDED USER, CONTRA-INDICATIONS AND PRECAUTIONS, RESIDUAL RISKS

- Essential performances: the acquisition/recording of the ECG, data transmission from the device to the App, and the saving of the data. During EMC disturbances the device can lose communication or stop acquisition for some seconds but restarts independently to work properly.
- The D-Heart device must be used on people of 18 years of age and above
- The acquisition of the ECG trace is done in autonomy without the intervention of a professional. The user, wherever they are and in total autonomy, can seek the assistance/opinion of a physician by providing an independently acquired ECG trace.

- Acquired traces must be evaluated by persons legally authorised by national laws and regulations, with appropriate skills in ECG analysis.
- The device is intended for supporting or providing useful information regarding the process of diagnosis or care of users at risk for or with heart diseases.
- The device is intended to be operated in hospital, in the general practitioner's office, in pharmacies and home environments.
- The conductive parts of electrodes must not come into contact with other conductive parts, including earth.
- Use ONLY disposable electrodes with REF SF01 from the manufacturer Shanghai INTCO Electrode Manufacturing Co.,Ltd or REF ES GS 3648 from the manufacturer Esvico Sp.z o.o. Contact info@d-heartcare.com to be informed on how to purchase the disposable electrodes.
- The device shall not be operated in combination with a cardiac defibrillator.
- The device is classified as Class B according to CISPR 11:2009 (Radio-frequency disturbance characteristics).
- The device is not suitable for direct cardiac application.
- An automatic check made by the D-Heart Portable ECG Device App advises and prevents the operator from situations when the device is inoperable (e.g. not all the electrodes are connected properly to the patient's body). This option is only available on iOS devices.
- The minimum operating time is 24 hours provided that the battery is new and fully charged.
- The battery charge time from depletion to 90 % charge in normal use and battery conditions is 2 hours.
- In order to charge the battery, the user must place the device inside its case and on the wireless charger included in the package. Battery charging is made by wireless technology.
- Level of battery charge can be checked by pushing the green icon on the top left side of the home page
- Don't connect the device to the body while the battery is being charged (please note that the device is not operative while the battery is being charged).
- After the charging has completed, leave the device getting cool before putting it in contact with the patient's body.
- The battery charging level must be checked by the user before operating the device by using the proper indication in the App.

- The isoelectric segments within the QRS complex are excluded from the Q-, R- or S-waves.
- The device is built in with specific filter settings to remove low-frequency components such as motion artefact, respiratory variation, and baseline wander.
- Before executing any measurement, check the device to make sure that there is no visible damage that may affect user's safety and measurement performance. Stop using the unit when there is obvious damage.
- The device only weighs 196gr and falling from the patient's hands will not harm them
- For D-Heart equipment that includes RF transmitters, the technical description shall include each frequency or frequency band of transmission, the type and frequency characteristics of the modulation and the EFFECTIVE RADIATED POWER.
 - Bluetooth low energy 4.1
 - Tx power: +4dB
 - Rx sensitivity: -88dB
 - Operating supply voltage: 1.7V – 3.6V
 - Operating temperature range: -40°C - +85°C
 - FHSS function
 - Modulation: GFSK
 - RF frequency: 2402MHz – 2480MHz
 - Operating channel bandwidth: 2MHz

14. MAINTENANCE, CLEANING AND TROUBLESHOOTING

a. CLEANING AND STERILIZING

- Turn off the device before cleaning. Wipe the device with a dry and clean cloth for cleaning. Do not allow any liquid to get inside the device.
- NEVER re-use the disposable electrodes.

b. MAINTENANCE

- Do not open the device case to avoid damages to internal components.
- Prevent any liquid from getting in the device as it would affect the safety and performance of the device.
- Software and firmware upgrade are notified to the user and automatically performed via the Smartphone App.
- The D-Heart device does not need any calibration
- The life time of the device is estimated to be 250 recharging cycles, which are estimated to equal to 4 years of usage.

c. TROUBLESHOOTING

c.1 TROUBLE: The device doesn't switch on

Possible reason	Solution
The battery is drained or almost drained.	Charge the battery.
The device is damaged.	Consult the following web site: https://www.d-heartcare.com/contact .

c.2 TROUBLE: The device doesn't connect to the Smartphone

Possible reason	Solution
The Bluetooth function of the Smartphone is switched off.	Switch the Bluetooth function of the Smartphone ON.
The Smartphone is connected to another D-Heart device.	Switch the other D-Heart OFF.
The device is damaged.	Consult the following web site: https://www.d-heartcare.com/contact .

c.3 TROUBLE: The automatic check in the App shows one or more disconnected electrodes

Possible reason	Solution
Disposable electrodes don't contact body well.	Place the disposable electrodes correctly.
The plugs are not connected to the disposable electrodes.	Connect the plugs to the disposable electrodes.
The device is damaged.	Consult the following web site: https://www.d-heartcare.com/contact .

c.4 TROUBLE: the device fails measuring the heart rate or the trace shows strong irrelevant waveforms

Possible reason	Solution
Disposable electrodes don't not contact body well.	Place the disposable electrodes correctly.
The plugs are not connected to the disposable electrodes.	Connect the plugs to the disposable electrodes.
Movement when measuring	When measuring, please keep quiet and avoid moving
Electromagnetic interference.	Keep away from interference source.
The device is damaged.	Consult the following web site: https://www.d-heartcare.com/contact .

c.5 TROUBLE: the battery doesn't charge

Possible reason	Solution
Wrong position on the wireless charger.	Place the device correctly on the wireless charger (the LED blinks slowly).
The device is damaged.	Consult the following web site: https://www.d-heartcare.com/contact .

15. LEGENDA OF SYMBOLS

The D-Heart device is certified according to the relevant regulations set by the European Union for electromedical devices (MDR (UE) 2017/745). The number "1370" identifies the notified body that verifies the compliance of the device to the applicable essential requirements.

Separate collection for Waste of Electric and Electronic Equipment (WEEE)

The D-Heart Portable ECG Device supports the Bluetooth® technology

The D-Heart Portable ECG Device is protected against water and dust with the level of protection IP22 as defined by the IEC standard 60529:

- The enclosure provides protection to access to hazardous parts for Fingers or similar objects
- The enclosure provides protection to Dripping water when tilted at 15°

Manufacturer

The D-Heart Portable ECG Device is classified as TYPE CF regarding the protection against electrical shock according to standard CEI EN 60601-2-25

Temperature limit

Humidity limitation

Atmospheric pressure limitation

Keep dry: store and use the device in a dry environment

“MR-unsafe”: Do not expose the device to a magnetic resonance (MR) environment

Medical Device

Date of Manufacture

Catalogue Number

Data Matrix (Example) describing concatenated UDI-DI and UDI-PI

Serial Number of each device

Unique Device Identifier: Concatenated UDI-DI and UDI-PI

Caution: it is mandatory to read the instructions for use

Instructions for use provided in paper format.
Reading the instructions for use is mandatory.

D-Heart App version

Internal Battery Electrical Potential

Model Number

Triman symbol

Non-corrugated cardboard packaging - Dispose of according to local laws and regulations

16. SUPPORT IN CASE OF NEED AND ADVERSE EVENT REPORTING

Please consult the following web site: <https://www.d-heartcare.com/contact> for assistance, if needed, in setting up, using or maintaining the D-Heart Portable ECG Device or to report unexpected operations or events. Incidents shall be reported. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority. Reporting form and information (dedicated exclusively to healthcare professionals) can be found at https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

Incidents related to D-Heart Portable ECG device shall also be reported to Strena Medical S.r.l on +39 010 3017000 or support@d-heartcare.com

17. ELECTROMAGNETIC COMPATIBILITY

This device is classified in class B according to IEC60601-1-2. This instrument has been tested and found to comply with the limits for medical devices to the IEC60601-1-2. These limits are intended to provide reasonable protection against harmful interference in a typical medical environment. This instrument generates, uses and can radiate radio frequency energies and, if not put in service and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular location. If the system does cause harmful interference to other devices, which can be determined by turning the system off and on, try to eliminate the interference by adopting one or more of the following measures:

- reorient and/or relocate the receiving device;
- increase the distance between the devices;
- connect the system to an outlet on a different circuit than that to which the other devices is connected;
- consult the manufacturer or field service technician for help.

Essential performance: continuous operation (measurement status).

a. SPECIFICATIONS AND TECHNICAL INFORMATION: EMC (ELECTROMAGNETIC COMPATIBILITY)

Guidance and manufacturer's declaration - electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

Guidance and manufacturer's declaration - electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity test	Compliance level
ETSI EN 301 489-1 V2.2.3 (2019-11) Electromagnetic compatibility and Radio spectrum Matters (ERM) Electromagnetic Compatibility (EMC) standard for radio equipment and services Part 1: Common technical requirements	10 V/m 80 MHz to 2.7 GHz
ETSI EN 301 489-17 V3.2.4 (2020-09) Electromagnetic compatibility and Radio spectrum Matters (ERM) Electromagnetic Compatibility (EMC) standard for radio equipment Part 17: Specific conditions for Broadband Data Transmission Systems	
CEI EN 60601-1-2:2016-04 EN 60601-1-2:2015-09 IEC 60601-1-2:2014-02 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances - Requirements and tests	
ETSI EN 301 489-1 V2.2.3 (2019-11) ElectroMagnetic Compatibility (EMC) standard for radio equipment and services; Part 1: Common technical requirements	3 V/m 2.7 GHz to 6.0 GHz
ETSI EN 301 489-17 V3.2.4 (2020-09) ElectroMagnetic Compatibility (EMC) standard for radio equipment and services; Part 17: Specific conditions for Broadband Data Transmission Systems	
CEI EN 60601-1-2:2016-04 EN 60601-1-2:2015-09 IEC 60601-1-2:2014-02 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances - Requirements and tests	

Guidance and manufacturer's declaration - electromagnetic immunity

Immunity test	Compliance level
Electrostatic Discharge	According to IEC 61000-4-2 ± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Rated Power magnetic fields	According to IEC 61000-4-8 30 A/m 50 Hz and 60 Hz
Proximity fields from RF wireless communications equipment	According to IEC 61400-4-3 CEI EN 60601-1-2:2016-04 Table 9 test level applied
Magnetic proximity fields IEC 61000-4-39	30 KHz, CW, 8 A/m 134.2 KHz, PM, 65 A/m 13560 KHz, PM, 7.5 A/m

Guidance and manufacturer's declaration – electromagnetic emissions

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Emission test	Compliance criteria
ETSI EN 300 328 v2.2.2 (2019-07) Data transmission equipment operating in the 2,4 GHz ISM band and using wide band modulation techniques; Harmonised Standard covering the essential requirements of article 3.2 of the Directive 2014/53/EU	All requirements are met

18. DISPOSAL OF THE DEVICE

According to Directive 2012/19/EU of the European Parliament and of the Council, consolidated text of 04/07/2018 on waste electrical and electrical equipment (WEEE), and Italian legislative decree 49/2014.

The crossed-out wheeled bin symbol printed either on the equipment or on the package indicates that the product shall be disposed separately from other waste at the end of its life cycle and that the user shall deliver the old equipment to collection centers authorized for electrical waste. Alternatively, the old equipment shall be sent back at the time of the purchasing of a new equipment of the same model in one-to-one proportion.

The correct disposal of waste according to the instructions above helps in preventing potential drawbacks for the environment and the public health and helps the recycling of the materials the equipment is made of. The incorrect disposal of the equipment is subject to administrative sanctions according to the existing law.

For further information do not hesitate to contact us at the details provided below.

19. CONTACT US

If you are in need of support please write to us via our web chat on our website: <https://www.d-heartcare.com/> write us an email at support@d-heartcare.com or call us at +39 010 30 17 000.

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VAT Number IT02335950990

www.d-heartcare.com

The Bluetooth® word mark and logos are registered trademarks of Bluetooth SIG, Inc. and any use of such marks by Strena Medical s.r.l. is under license. Google Play and the Google Play logo are trademarks of Google Inc. Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

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Copyright 2024 Strena Medical S.r.l. SB
Strena Medical S.r.l. Società Benefit
Strena Medical S.r.l. SB