ri-mini[®] pen-scope[®]

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1. Important information to take note of before taking the product into operation

You have acquired a valuable **Riester** diagnostic set manufactured in compliance with Directive 93/42/EEC for medical products and subject to continuous stringent quality control, whose excellent quality will ensure reliable diagnoses. Please read these Operating Instructions carefully prior to startup and keep in a safe place. Should you have any queries, please contact the Company or your **Riester** Agent who will be pleased to assist you. For addresses see last page of these Operating Instructions. The address of your authorised **Riester** Agent will be supplied to you on request. Please note that any instruments described in these Operating Instructions are only suited for application by trained operators. Please also note that correct and safe operation of instruments will only be guaranteed when **Riester** instruments and accessories are used throughout.

Warning:

Please note that the proper and safe functioning of our instruments is only guaranteed if both the instruments and their accessories are exclusively from Riester. The use of other accessories may result in increased electromagnetic emissions or reduced electromagnetic immunity of the device and may lead to incorrect operation.

Caution / Contraindications

- There may be a risk of ignition of gases when the instrument is used in the presence of flammable mixtures or mixtures of pharmaceuticals.
- The instrument heads and battery handles must never be placed in liquids.
- In the case of a longer eye examination using the ophthalmoscope, the intensive exposure to light can damage the retina.
- The product and the ear specula are non-sterile. Do not use on injured tissue.
- Use new or sanitized specula to limit the risk of cross-contamination.
- The disposal of used ear specula must occur in accordance with current medical practices or local regulations regarding the disposal of infectious, biological medical waste.
- Only use accessories / consumables from Riester or Riester-approved accessories / consumables.
- Cleaning frequency and sequence must comply with regulations for the cleaning of non-sterile products in their respective facility. Cleaning / disinfection instructions in the operating manual must be observed.
- The product may only be used by trained personnel.

Safety instructions:

Manufacturer

CE marking

**Temperature limits
in °C for storage and transport**

**Temperature limits
in °F for storage and transport**

Relative humidity

Fragile, handle with care

Store in a dry place

„Green Dot“ (country-specific)

Warning, this symbol indicates a potentially dangerous situation

Device of protection class II

Application part type B

For single use only

Attention: Used electrical and electronic equipment should not be treated as normal household waste but should be disposed of separately in accordance with national or EU directives

Batch code

Serial number

Please observe the operating instructions

2.1. Purpose

Riester battery handles described in these Instructions for Use supply the instrument heads with power (the lamps are included in appropriate instrument heads), also serving as a bracket.

2.2. Readiness for operation (insertion and removal of batteries)

Turn off instrument head from handle in counter-clockwise direction. Insert two commercial type "AA" Mignon alkaline batteries of 1.5 V (IEC standard reference LR6) into the case of the handle with the plus poles towards the upper section of the handle.

Warning:

- Should the unit not be used for an extended period of time or whilst travelling, remove batteries from handle.
- Insert new batteries when light intensity of the unit is reduced, thus affecting examination.
- For maximum light yield it is recommended to always insert new high-quality batteries on replacement.
- Ensure that no fluid or condensation penetrates into the handle. Disposal Please note that batteries are subject to separate disposal. For details ask your local authority and/or your environmental officer.

Disposal

Please note that batteries are subject to separate disposal. For details ask your local authority and/or your environmental officer.

2.3. Attachment of instrument heads

Turn instrument head in clockwise direction on to the handle. The instrument head may be fixed in the required position by turning the ring on the instrument head up to the stop.

2.4. Starting and stopping

The clip on the handle is equipped with an On/Off slide. When pushing the slide up, the unit is switched on, when pushing it down, the unit is off.

2.5 Instructions for care

General information

Cleaning and disinfection of the medical devices serves to protect the patient, the user and third parties and to preserve the value of the medical devices. Due to the product design and the materials used, no defined limit can be specified for the maximum number of reprocessing cycles that can be carried out. The life span of the medical devices is determined by their function and by gentle handling of the devices. Defective products must undergo the reprocessing procedure described before being returned for repair.

Cleaning and disinfection

The battery handles can be cleaned externally with a moist cloth until visually clean. Wipe-disinfection as specified by the disinfectant manufacturer. Only disinfectants with proven efficacy should be used, taking into account the national requirements. After disinfection, wipe the instrument down with a moist cloth to remove possible disinfectant residues.

PLEASE NOTE!

Never immerse the handles in liquids! Take care to ensure that no liquids get inside the casing! This item is not approved for automated reprocessing and sterilization. These procedures cause irreparable damage!

3. Otoscope and accessories

3.1. Purpose

Riester otoscopes described in these Instructions for Use have been produced for lighting and examination of the auditory canal, combined with a **Riester** ear speculum.

3.2. Insertion and removal of ear speculum
Position the selected speculum on the chromium plated metal socket of the otoscope. Turn speculum to the right until a resistance is felt. The size of the speculum is marked on the reverse.

3.3. Introduction of external instruments into the ear
When intending to introduce external instruments into the ear (such as forceps), turn magnifying glass (approx. 2.5 X enlargement) on otoscope head in counterclockwise direction. Replace cover lens in reverse direction.

3.4. Pneumatic test
In order to perform a pneumatic test (= examination of the ear drum), you will require a bulb which is not included in the normal scope of supply but may be ordered separately (see Spare parts and accessories). Take metal connector which is not included in the normal scope of supply but may be ordered separately (see Spare parts and accessories) and insert in recess provided on the side of the otoscope head. Attach hose of bulb to connector. Carefully introduce the required air volume into the auditory canal.

3.5. Replacement of lamp
ri-mini® Remove instrument head from battery handle. The lamp is in the bottom section of the instrument head. Remove lamp, by using your thumb and forefinger or a suitable tool, from instrument head. Firmly insert new lamp.

pen-scope®
Screw off the metal fitting. Release lamp by turning in counter-clockwise direction. Insert new lamp by screwing in in clockwise direction and screw on the metal fitting again.

3.6. Technical data of the lamp			
ri-mini® otoscope:	XL, 2.5 V	750 mA	mean life span 16.5h
pen-scope® otoscope:	vacuum, 2.7 V	300 mA	mean life span 15h
	XL, 2.5 V	750 mA	mean life span 16.5h

3.7 Instructions for care
General information
Cleaning and disinfection of the medical devices serves to protect the patient, the user and third parties and to preserve the value of the medical devices. Due to the product design and the materials used, no defined limit can be specified for the maximum number of reprocessing cycles that can be carried out. The life span of the medical devices is determined by their function and by gentle handling of the devices. Defective products must undergo the reprocessing procedure described before being returned for repair.

Cleaning and disinfection
The otoscope can be cleaned externally with a moist cloth until visually clean. Wipe-disinfection as specified by the disinfectant manufacturer. Only disinfectants with proven efficacy should be used, taking into account the national requirements. After disinfection, wipe the instrument down with a moist cloth to remove possible disinfectant residues.

PLEASE NOTE! 
Never immerse the otoscope in liquids! Take care to ensure that no liquids get inside the casing! This item is not approved for automated reprocessing and sterilization. These procedures cause irreparable damage!

Sterilization
a) Reusable ear specula
The ear specula can be sterilized in the steam sterilizer at 134°C with 10 minutes hold time.

b) Single use ear specula
 For single use only

ATTENTION: Repeated use could cause infection

3.8. Spare parts and accessories

Reusable ear specula

Item no. 10775	Pack of 10 pcs	2 mm, black
Item no. 10779	Pack of 10 pcs	2,5 mm, black
Item no. 10783	Pack of 10 pcs	3 mm, black
Item no. 10789	Pack of 10 pcs	4 mm, black
Item no. 10795	Pack of 10 pcs	5 mm, black

Single use ear specula

Item no. 14061-532	Pack of 100 pcs	2 mm, black
Item no. 14062-532	Pack of 500 pcs	2 mm, black
Item no. 14063-532	Pack of 1000 pcs	2 mm, black
Item no. 14061-531	Pack of 100 pcs	2,5 mm, black
Item no. 14062-531	Pack of 500 pcs	2,5 mm, black
Item no. 14063-531	Pack of 1000 pcs	2,5 mm, black
Item no. 14061-533	Pack of 100 pcs	3 mm, black
Item no. 14062-533	Pack of 500 pcs	3 mm, black
Item no. 14063-533	Pack of 1000 pcs	3 mm, black
Item no. 14061-534	Pack of 100 pcs	4 mm, black
Item no. 14062-534	Pack of 500 pcs	4 mm, black
Item no. 14063-534	Pack of 1000 pcs	4 mm, black
Item no. 14061-535	Pack of 100 pcs	5 mm, black
Item no. 14062-535	Pack of 500 pcs	5 mm, black
Item no. 14063-535	Pack of 1000 pcs	5 mm, black

Replacement lamps

For pen-scope® otoscope

No. 10488	Pack with 6 vacuum lamps, 2.7 V
No. 10489	Pack with 6 xenon lamps, 2.5 V

For ri-mini® otoscope

Item No. 10600	Pack with 6 xenon lamps, 2.5 V
----------------	--------------------------------

Other spare parts

No. 10960	Bulb for pneumatic test
No. 10961	Connector for pneumatic test
No. 10963	Magnifying glass for ri-mini® and pen-scope® Otoscope, black

4. Ophthalmoscope and accessories

4.1. Purpose

Riester ophthalmoscopes described in these Instructions for Use have been designed for the examination of the eye and its background.

ATTENTION!

Because prolonged intense exposure to light can damage the retina, the use of the eye exam device should not be unnecessarily prolonged, and the brightness setting should not be set higher than needed for a clear representation of the target structures.

The irradiation dose of the photochemical exposure to the retina is the product of irradiance and duration of irradiation. If the irradiance is reduced by half, the irradiation time may be twice as long to reach the maximum limit.

Although no acute optical radiation hazards have been identified for direct or indirect ophthalmoscopes, it is recommended that the intensity of light directed into the patient's eye be reduced to the minimum required for examination/diagnosis. Infants / children, aphasics and people with eye diseases are at a higher risk. The risk may be increased if the patient has already been examined with this or another ophthalmological instrument during the last 24 hours. This is especially true when the eye has been exposed to retinal photography.

The light of this instrument may be harmful. The risk of eye damage increases with the duration of irradiation. An irradiation period with this instrument at maximum intensity of longer than →5 min. exceeds the guideline value for hazards.

This instrument does not pose a photobiological hazard according to DIN EN 62471 but still features a safety shutdown after 2 / 3 minutes.

4.2. Lens wheel and correcting lenses

The correcting lenses may be adjusted on the lens wheel. The following correcting lenses are available:

diopeters 0 to + 20 and 0 to - 20. Readings will be displayed on a lit panel. Plus values are displayed in black digits, minus values in red digits.

4.3. Apertures and filters

The following apertures and/or filters may be selected by the aperture and filter wheel:

Aperture	Function
Semi-circle:	for examinations with turbid lenses
Small circle:	for reduction of reflexes of small pupils
Large circle:	for standard fundus examination
Fixation star:	for definition of central and eccentric fixation

The following filters are only available in the **ri-mini®** ophthalmoscope: Red-free filter To increase contrast for assessment of (green filter) changes in fine vessels, i.e. retinal haemorrhages

4.4. Lamp replacement

Remove instrument head from battery handle. The lamp is in the bottom section of the instrument head. Remove lamp, by using your thumb and forefinger or a suitable tool, from the instrument head. Insert new lamp with the pin on the lamp fitted in the recess provided on the head of the case.

4.5 Technical data of the lamp

ri-mini®	Ophthalmoscope XL, 2.4 V	750 mA	mean life span 15h
pen-scope®	Ophthalmoscope XL, 2.5 V	750 mA	mean life span 15h

4.6 Instructions for care

General information

Cleaning and disinfection of the medical devices serves to protect the patient, the user and third parties and to preserve the value of the medical devices. Due to the product design and the materials used, no defined limit can be specified for the maximum number of reprocessing cycles that can be carried out. The life span of the medical devices is determined by their function and by gentle handling of the devices. Defective products must undergo the reprocessing procedure described before being returned for repair.

Cleaning and disinfection

The ophthalmoscope can be cleaned externally with a moist cloth until visually clean. Wipe-disinfection as specified by the disinfectant manufacturer. Only disinfectants with proven efficacy should be used, taking into account the national requirements. After disinfection, wipe the instrument down with a moist cloth to remove possible disinfectant residues.

PLEASE NOTE!

Never immerse the ophthalmoscope in liquids! Take care to ensure that no liquids get inside the casing! This item is not approved for automated reprocessing and sterilization. These procedures cause irreparable damage!

4.7 Spare parts

for **pen-scope®** ophthalmoscope

No. 10473 Pack with 6 xenon lamps, 2.5 V

for **ri-mini®** ophthalmoscope

No. 10605 Pack with 6 xenon lamps, 2.5 V

5. Maintenance

These instruments and their accessories do not require any specific maintenance. Should an instrument have to be examined for any specific reason whatsoever, please return it to the Company or an authorised **Riester** dealer in your area. Addresses to be supplied on request.

6. Notices

Ambient temperature:	0 ° to +40 ° C
Relative Humidity:	30% to 70% noncondensing
Storage location:	-10° to +55°
Relative Humidity:	10% to 95%

7. Electromagnetic compatibility

Accompanying documents according to IEC 60601-1-2, 2014, Ed. 4.0

Attention:

Medical electrical equipment is subject to special precautions regarding electromagnetic compatibility (EMC).

Portable and mobile radio frequency communication devices can affect medical electrical equipment. The ME device is for operation in an electromagnetic environment for home health care and intended for professional facilities such as industrial areas and hospitals.
The user of the device should ensure that it is operated within such an environment.

Warning:

The ME device may not be stacked, arranged or used directly next to or with other devices. When operation is required to be close to or stacked with other devices, the ME device and the other ME devices must be observed in order to ensure proper operation within this arrangement. This ME device is intended for use by medical professionals only. This device may cause radio interference or interfere with the operation of nearby devices. It may become necessary to take appropriate corrective measures, such as redirecting or rearranging the ME device or shield.
The rated ME device does not exhibit any basic performance features in the sense of EN60601-1, which would present an unacceptable risk to patients, operators or third parties should the power supply fail or malfunction.

Warning:

Portable RF communications equipment (radios), including their accessories, such as antenna cables and external antennas, should not be used less than 30 cm (12 inches) away from the ri-mini® and pen-scope® instrument head parts and cables specified by the manufacturer. Failure to comply may result in a reduction in the device's performance features.

Guidance and manufacturer's declaration - electromagnetic emission		
The ri-mini® and pen-scope® instruments are intended for use in the electromagnetic environment specified below. The customer or user of the ri-mini® and pen-scope® should ensure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment – guidance
RF emissions HF- emissions pursuant to CISPR 11	Group 1	The ri-mini® and pen-scope® use RF energy exclusively for internal functions. Therefore, its RF emissions are very low and unlikely to disturb nearby electronic devices. The ri-mini® and pen-scope® are intended for use in all establishments, including residential areas and those directly connected to a public supply network that also supplies buildings used for residential purposes.
RF emissions HF- emissions pursuant to CISPR 11	Class B	
Emissions of harmonics IEC 61000-3-2	Not applicable	
Emissions of voltage fluctuations, flicker IEC 61000-3-3	Not applicable	

Directives and manufacturer's declaration - Electromagnetic immunity			
The ri-mini® and pen-scope® instruments are intended for use in the electromagnetic environment specified below. The customer or user of the ri-mini® and pen-scope® should ensure that it is used in such an environment.			
Immunity testing	IEC 60601 test level	Compliance	Electromagnetic environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	Con: ±8 kV Air: ±2,4,8,15 kV	Con: ±8 kV Air: ±2,4,8,15 kV	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Fast transient electrical disturbances / bursts IEC 61000-4-4	5/50 ns, 100 kHz , ±2 kV	Not applicable	The quality of the supply voltage should be that of a typical business or hospital environment.
Surge voltage IEC 61000-4-5	± 0.5 kV voltage Phase-to-phase conductor ± 2 kV voltage Line-to-earth	Not applicable	The quality of the supply voltage should be that of a typical business or hospital environment.
Voltage dips, short-term interruptions and fluctuations in the supply voltage acc. to IEC 61000-4-11	<0% UT 0.5 period at 0.45, 90, 135, 180, 225, 270 and 315 degrees 0% UT 1 period and 70% UT 25/30 periods single-phase: at 0 degrees (50/60 Hz)	Not applicable	The quality of the supply voltage should be that of a typical business or hospital environment.
Magnetic field with efficiency-rated frequencies IEC 61000-4-8	30A/m 50/60 Hz	30A/m 50/60 Hz	Mains frequency magnetic fields should be at a level characteristic of a typical location in a typical commercial hospital environment.
NOTE UT is the AC source. Mains voltage before the application of the test level.			

Directives and manufacturer's declaration - Electromagnetic immunity			
The ri-mini® and pen-scope® instruments are intended for use in the electromagnetic environment specified below. The customer or user of the ri-mini® and pen-scope® should ensure that it is used in such an environment.			
Immunity testing	IEC 60601 test level	Compliance	Electromagnetic environment - Instructions
Guided RF Disturbances acc. to IEC61000-4-6	3 Vrms 0.5 MHz bis 80MHz 6 V in ISM frequency bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz	Not applicable	Portable and mobile RF communications equipment should not be used closer to any part of the non-contact ri-mini® and pen-scope®, including the cables, than the recommended distance, which is calculated using the equation applicable to the transmitter frequency. Recommended separation distance $d = 1.2 \times P$ 80 MHz to 800 MHz $d = 2.3 \times P$ 800 MHz to 2.7 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b
			Interference may occur in the vicinity of equipment marked with the following symbol:
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	10 V/m	
Proximity fields from RF wireless communications equipment	380 - 390 MHz 27 V/m; PM 50%; 18 Hz 430 - 470 MHz 28 V/m; (FM ±5 kHz, 1 kHz sine) PM; 18 Hz 704 - 787 MHz 9 V/m; PM 50%; 217 Hz 800 - 960 MHz 28 V/m; PM 50%; 18 Hz 1700 - 1990 MHz 28 V/m; PM 50%; 217 Hz 2400 - 2570 MHz 28 V/m; PM 50%; 217 Hz 5100 - 5800 MHz 9 V/m; PM 50%; 217 Hz	27 V/m 28 V/m 9 V/m 28 V/m 28 V/m 9 V/m	
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a) field strengths of fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM broadcast and television transmission cannot be theoretically predicted with accuracy. In order to assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic survey should be considered. If the measured field strength at the location where the ri-mini® and pen-scope® are used exceeds the above-mentioned RF degree of compliance, the ri-mini® and pen-scope® should be observed to ensure normal operation. If abnormal performance is observed, additional measures may be required, such as reorienting or moving the ri-mini® and pen-scope®.			
b) With a frequency range over 150 kHz to 80 MHz, the field strengths should be less than 3 V/m.			

Recommended distances between portable and mobile RF communications equipment and the ri-mini® and pen-scope®

The ri-mini® and pen-scope® are intended for use in an electromagnetic environment in which RF emissions are controlled. The customer or user of the ri-mini® and pen-scope® can help to avoid electromagnetic interference by observing the minimum distance between portable and mobile RF communications equipment (transmitters) and the ri-mini® and pen-scope® in accordance with the maximum output power of the communication equipment.

Rated maximum output power of the transmitter (W)	Separation distance according to the frequency of the transmitter (m)			800 MHz to 2.7 GHz
	150 KHz to 80 MHz		80 MHz to 800 MHz	
0.01	0.12	0.12	0.23	0.23
0.1	0.38	0.38	0.73	0.73
1	1.2	1.2	2.3	2.3
10	3.8	3.8	7.3	7.3
100	12	12	23	23
For transmitters with a maximum output power not listed above, the recommended distance d in metres (m) can be estimated using the equation for the transmitter frequency, where P is the maximum output power of the transmitter in watts (W), according to the transmitter manufacturer.				
NOTE 1 At 80 MHz and 800 MHz, the separation distance applies to the higher frequency range.				
NOTE 2 These guidelines may not apply in all situations. The electromagnetic propagation is affected by absorption and reflection of structures, objects and people.				

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WARRANTY

This product has been manufactured under the strictest quality standards and has undergone a thorough final quality check before leaving our factory.

We are therefore pleased to be able to provide a warranty of **2 years from the date of purchase** on all defects, which can verifiably be shown to be due to material or manufacturing faults. A warranty claim does not apply in the case of improper handling.

All defective parts of the product will be replaced or repaired free of charge within the warranty period. This does not apply to wearing parts.

For **r1 shock-proof**, we grant an additional warranty of 5 years for the calibration, which is required by CE-certification.

A warranty claim can only be granted if this Warranty Card has been completed and stamped by the dealer and is enclosed with the product.

Please remember that all warranty claims have to be made during the warranty period. We will, of course, be pleased to carry out checks or repairs after expiry of the warranty period at a charge. You are also welcome to request a provisional cost estimate from us free of charge.

In case of a warranty claim or repair, please return the **Riester** product with the completed Warranty Card to the following address:

Rudolf Riester GmbH
Dept. Repairs RR
Bruckstr. 31
72417 Jungingen
Germany

Serial number or batch number Date, Stamp and signature of the specialist dealer

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